

## Choice<sup>®</sup> Cleanroom Socks

99% Continuous Filament Polyester / 1% Spandex Cleanroom Socks

Choice<sup>®</sup> Cleanroom Socks are made of continuous filament polyester yarn and are laundered and packaged in a Class 5 cleanroom to minimize bioburden levels and reduce particulate contamination in critical environments. The mid-calf length and elastic cuff result in a comfortable secure fit for extended wear.



### Features

- 99% Continuous filament polyester, 1% spandex
- One size fits all
- Mid-calf length with elastic cuff for secure fit
- Soft and comfortable
- Laundered and packaged in an ISO Class 5 cleanroom
- Color: White
- 10 pairs per pack

### Technology

- Constructed with continuous filament yarns to reduce fiber and particle shedding
- Cleanroom laundered and packaged

### Advantages

- Meets the requirements of EU GMP Annex 1 section 7.15
- Cleanroom laundry process reduces fiber and particle shedding and minimizes bioburden
- Over 1000 times cleaner than regular polyester/cotton socks
- Comfortable and durable for all day wear
- A cost effective, clean and disposable alternative to operators' own socks
- Autoclavable (121°C, 30 minutes)
- Bulk packaging available

### Technical Data

Description	Test Method	Average Data
Polyester/cotton sport sock	Helmke Drum part $\geq 0.5\mu\text{m} / \text{ft}^3 / \text{sock}$	45,000
Choice <sup>®</sup> Cleanroom Sock	Helmke Drum part $\geq 0.5\mu\text{m} / \text{ft}^3 / \text{sock}$	120

**Note:** Data shown are typical average values

### Order Information:

Product	Number	Prs/Pk	Pks/Cs	Style Packed
Choice <sup>®</sup> Cleanroom Sock	CHSOCKR1050	10	50	Stacked

[www.berkshire.com](http://www.berkshire.com)

Contact: Tel 1 800 242 7000 / 1 413 528 2602  
info@berkshire.com

America	Tel 1 413 528 2602	info@berkshire.com
Europe	Tel + 44 1953 562800	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	enquiries@berkshire.com.sg
Japan	Tel 81 3 4530 9883	master@berkshire.co.jp